

DIAGNOSTIC CONSUMABLES

Diagnostic Reagent Sample Kits

As diagnostic tests move closer to the point-of-care, reagent storage and delivery methods become a challenge. J-Pac Medical's lab-on-chip reagent blisters allow reagents to be stored directly on the microfluidic cartridge, enabling the smaller footprint required by point-of-care diagnostic tests.

Contact customer service at 1-603-692-9955 to order sample kits.

About J-Pac Medical

J-Pac Medical is the trusted manufacturing and packaging outsourcing partner to medical device companies seeking to deliver superior quality, improve time-to-market, and simplify the supply chain for single-use medical devices.

Dimensions	Sample Blister Fluid Volume			
	30 μ l	130 μ l	500 μ l	1,000 μ l
Height	0.108	0.175	0.26	0.268
Diameter	0.394	0.585	0.875	0.988
Flange	0.065	0.118	0.096	0.118
Pierce Dia.	0.264	0.375	0.625	0.752
Quantity	Water Filled (In-Stock)			
25	\$250			
100	\$500			
250	\$1,000			
500	\$1,500			
Quantity	Customer Supplied Reagent			
<1000	\$1500 plus \$5/ea			
1000-2000	\$1500 plus \$4/ea			

